



## Infinite® F50 – Reliable. Validated. Compact.

New state-of-the-art 8-channel absorbance reader with latest LED technology

---

## Infinite F50 – Innovative design with top price-performance ratio

The Infinite F50 is Tecan's next-generation, state-of-the-art 8-channel absorbance microplate reader that provides not only accurate, reproducible and fast measurements but sets a new standard in compact, innovative design. With its easy-to-use control and data analysis software, the Infinite F50 is ideal for a variety of ELISA applications.

### Reliable and validated hardware and software in an outstandingly compact footprint

#### Maintenance-free

The Infinite F50 uses the latest LED technology to replace the traditional halogen lamp light source, providing longer maintenance-free life, lower energy consumption and compact dimensions. LEDs typically last 10 times longer than halogen lamps at the same efficiency, making inconvenient lamp changes and waiting for the lamp to warm up a thing of the past.

#### Proven reliability

Reliability of the Infinite F50 is proven by Tecan's unique quality tests resulting in a Mean Time to First Failure (MTTF) value of >100,000 plates, typically corresponding to several years of successful operation.

#### Warranty

The 24-month warranty period for the Infinite F50 underlines the quality of Tecan's detection instruments that is demonstrated by the installed base of more than 20,000 Tecan microplate readers worldwide.

#### Validated

Both the Infinite F50 absorbance reader and Magellan™ software have been designed, tested, verified and validated according to ISO 13485:2003 quality system standards. Together with the optional Magellan Tracker software, the system meets the 98/79/EC IVD directive for in vitro diagnostic products and the FDA regulation 21 CFR Part 11, ensuring technical compliance and integrity of generated data.

#### Compact

The innovative, state-of-the-art engineering enables unprecedented compactness, reflected in a footprint of about four microplates. This reader can easily fit into any laboratory, saving valuable bench space.

## Comprehensive and powerful data analysis

The Infinite F50 comes complete with Magellan – Tecan's established reader control and data analysis software – allowing quick and easy measurements, comprehensive data reduction and excellent data presentation on cost-effective netbooks, as well as on laptops or standard PCs.

### Powerful

- full qualitative and quantitative EIA analysis
- all major curve fittings (eg. point-to-point, linear and non-linear regression, polynomial, logit-log, four- and five-parameter fits)
- comprehensive QC function (such as plate-to-plate QC for trending of standards and controls displayed in a Levy-Jennings graph)
- a range of data import and export options (Excel®, ASCII, ASTM)
- a range of kinetic data reduction features (such as mean and max slope, enzyme kinetics, etc.)\*\*

### Easy access

- wizard-based user interface guides you through the software step-by-step
- available in eight languages (Chinese, English, French, German, Italian, Portuguese, Russian and Spanish)
- example files for qualitative and quantitative ELISA assays and learning tutorials

### Compliant

Magellan Tracker software provides electronic signatures, records and multilevel user management to offer full functionality to meet FDA's Title 21 CFR Part 11 and 98/79/EC IVD-D, allowing diagnostic applications in clinical laboratories and other regulated environments. The Infinite F50 is also compliant with the Restriction of Hazardous Substances (RoHS) and the Waste Electrical and Electronic Equipment (WEEE) directives.

## Infinite F50 features and benefits

- 8-channel absorbance reader for fast measurements
- long-life LED for maintenance-free operation
- outstanding reliability proven by MTTF benchmark value
- endpoint, kinetic and multi-label measurements for a variety of applications
- shaking with variable time and speed
- Magellan software for powerful and validated data analysis
- learning tutorials, example files and eight different languages for getting started easily
- 98/79/EC IVD-D and FDA 21 CFR Part 11 for regulated environments

## Available options

The hand-held barcode scanner (LIS option) is a time-saving and convenient option that allows reading of barcodes from sample IDs or plate IDs with automatic transfer into Magellan software and data export to the appropriate Laboratory Information System (LIS).



**m magellan**  
data analysis software

## Infinite F50 – Typical performance values

General specifications	
Detection mode	absorbance
Light source	LED, self-calibrating
Photodetector	silicon photodiodes
Plate type	96-well microplates (including strip-well microplates)
Measurement modes	endpoint, kinetic and multi-label measurements
Measurement channels	8
Measurement range	0 – 4 OD
Resolution	0.0001 OD
Wavelength range	400 – 750 nm
Wavelength selection	up to 8 filters on wheel standard filters 405, 450, 492 and 620 nm***
Accuracy at 450 / 492 nm	
0.000 – 2.000 OD	≤ (0.5 % + 0.010 OD) typical
2.000 – 3.000 OD	≤ (1 % + 0.010 OD) typical
Precision / Reproducibility at 450 / 492 nm	
0.000 – 2.000 OD	≤ (0.5 % + 0.005 OD)
2.000 – 3.000 OD	≤ (1.0 % + 0.005 OD)
Linearity	
0.000 – 2.000 OD	≤ 1.0 %
2.000 – 3.000 OD	≤ 1.5 %
Shaking	linear shaking in 4 different modes
Computer interface	USB 1.1 / 2.0
Dimensions	height: 13.4 cm (5.28 in) width: 34.7 cm (13.66 in) depth: 18.9 cm (7.44 in)
Weight	< 3 kg
Power supply	External power supply: 100 – 240 V AC, 50 / 60 Hz, max. 1.2 A (auto sensing)

## Tecan – Who we are

Tecan is a leading global provider of laboratory instruments and solutions in biopharmaceuticals, forensics, clinical diagnostics, academic centers and life science industries, and specializes in the development and production of automation solutions, detection instruments such as microplate readers, microarray-related products and washers.

Founded in Switzerland in 1980, Tecan has manufacturing, research and development sites in both North America and Europe, and maintains a sales and service network in 52 countries. To date, Tecan has distributed approximately 20,000 microplate readers worldwide and is committed to continuous technological improvements and compliance to global quality standards.

**Austria** +43 62 46 89 33 **Belgium** +32 15 42 13 19 **China** +86 10 5869 5936 **Denmark** +45 70 23 44 50 **France** +33 4 72 76 04 80 **Germany** +49 79 51 94 170  
**Italy** +39 02 92 44 790 **Japan** +81 44 556 73 11 **Netherlands** +31 18 34 48 174 **Portugal** +351 21 000 82 16 **Singapore** +65 644 41 886 **Spain** +34 93 490 01 74  
**Sweden** +46 31 75 44 000 **Switzerland** +41 44 922 89 22 **UK** +44 118 9300 300 **USA** +1 919 361 5200 **Other countries** +41 44 922 8125

\*together with Magellan Tracker software

\*\*For further details please refer to Magellan software features ([www.tecan.com/magellan](http://www.tecan.com/magellan)).

\*\*\*additional filters available separately

Tecan Group Ltd. makes every effort to include accurate and up-to-date information within this publication; however it is possible that omissions or errors might have occurred. Tecan Group Ltd. cannot, therefore, make any representations or warranties, expressed or implied, as to the accuracy or completeness of the information provided in this publication. Changes in this publication can be made at any time without notice.

Tecan and Infinite are registered trademarks and Magellan is a trademark of Tecan Group Ltd. Männedorf, Switzerland.

Excel is a registered trademark of Microsoft Corporation or its subsidiaries in the US and/or other countries.

© 2010, Tecan Trading AG, Switzerland, all rights reserved. 11003BR100759

[www.tecan.com/InfiniteF50](http://www.tecan.com/InfiniteF50)

